

Optical Nonlinear Spectroscopy of Gold and Silicon Nanoparticles


Sean M. Anderson

Centro de Investigaciones en Óptica, A.C


January 30, 2012



When you combine. . .




with these. . .





Historic use of metallic nanoparticles.¹


¹Image: Raimond Spekking.



Nanoparticles have large surface to volume ratios.


Our Samples²

- Created via ion implantation
- Implanted on Suprasil 300 substrate
- 8×8 and 8×4 mm transparent samples



²Provided by Dr. Alejandro Reyes and Dr. Alicia Oliver of the IF-UNAM.

Au Sample



TEM scan.


Specifications

- Size: ~ 3 nm
- Density: 5.8×10^{16} nps/cm³
- Dose: 3.1×10^{16} atoms/cm²
- Implantation Energy: 2 MeV

Au²⁺ Sample


Specifications

- Size: ~ 3 nm
- Density: 9.8×10^{16} nps/cm³
- Dose: 5.0×10^{16} atoms/cm²
- Implantation Energy: 2 MeV



TEM scan.

Si Sample




Specifications

- Size: ~ 20 nm
- Dose: 2.5×10^{17} atoms/cm²
- Implantation Energy: 1.5 MeV

TEM scan.

Centrosymmetric Materials

A centrosymmetric material is a material that displays inversion symmetry, such that $p(x, y, z) \rightarrow p(-x, -y, -z)$.




- Many nonlinear materials are centrosymmetric
- Nanospheres are definitely centrosymmetric
- The material in these nanoparticles is centrosymmetric

Second Harmonic Generation (SHG)

Characteristics³

- Two photons of the same frequency combine
- Create one photon of double the frequency





³Image: Jon Chui

Sum-frequency Generation (SFG)

Characteristics

- Two photons of different frequency combine
- Create one photon that is the sum of both frequencies



Second-order Nonlinear Effects

Early work^{4 5} demonstrated that second-order processes

- Are dipole forbidden in the bulk of centrosymmetric materials
- Are related to $\chi^{(2)}$, the nonlinear susceptibility
- Have bigger dipolar (surface) than quadrupolar contributions

Second-order processes are well studied for flat surfaces, but what about round materials like nanospheres?

⁴J.A. Armstrong et al. *Physical Review*, 127(6):1918–1939, Sep 1962.

⁵N. Bloembergen et al. *Physical Review*, 128(2):606–622, Oct 1962.

Nonlinear Response of Nanoparticles

- Early work^{6 7} show that nanospheres have both dipolar and quadrupolar contributions
- A few years later, Mochán et al.⁸ determine the nonlinear polarization for an array of nanospheres.

⁶J.I. Dadap et al. *Physical Review Letters*, 83(20):4045–4048, 1999.

⁷V.L. Brudny et al. *Physical Review B*, 62(16):11152, 2000.

⁸W.L. Mochán et al. *Physical Review B*, 68(8):085318, 2003.

Theory

The dipole moment for a single nanosphere is

$$\mathbf{p}^{(2)} = \gamma^e \mathbf{E}^{\text{ex}}(0) \cdot \nabla \mathbf{E}^{\text{ex}}(0) + \gamma^m \mathbf{E}^{\text{ex}}(0) \times [\nabla \times \mathbf{E}^{\text{ex}}(0)]. \quad (1)$$

The quadrupole moment for a single nanoshpere is

$$\mathbf{Q}^{(2)} = \gamma^q \left(\mathbf{E}^{\text{ex}}(0) \mathbf{E}^{\text{ex}}(0) - \frac{1}{3} [\mathbf{E}^{\text{ex}}(0)]^2 \mathbf{1} \right), \quad (2)$$

where the γ 's are unknown nonlinear response functions. $\mathbf{p}^{(2)}$ is nonlocal because of the field derivative, while $\mathbf{Q}^{(2)}$ is local.

The total nonlinear polarization for all the nanospheres is

$$\mathbf{P}^{nl} = n_s \mathbf{p}^{(2)} - \frac{1}{6} \nabla \cdot n_s \mathbf{Q}^{(2)}, \quad (3)$$

where n_s is the nanocrystal volume density. We substitute equations (1) and (2) to obtain

$$\mathbf{P}^{(2)} = \Delta' (\mathbf{E} \cdot \nabla \mathbf{E}), \quad (4)$$

where $\Delta' \equiv n_s (\gamma^e - \gamma^m - \gamma^q/6)$ and represents a kind of bulk response function.

Summary

Nonlinear response depends on

- Nonlocal excitation of the electric dipole moment
- Local excitation of the electric quadrupole moment
- The strength of the incident beam and
- The form (plane wave, Gaussian beam, polarization, etc.)
- The quadrupolar $(\mathbf{E} \cdot \nabla) \mathbf{E}$ term

What's the best way to enhance this signal?

The XP2SHG/SFG Technique

- Early work⁹ shows that using two cross-polarized beams reduces the number of unknown $\chi^{(2)}$ components
- Wang et al. manage to discern surface and bulk contributions using two beams¹⁰
- Using two beams greatly increases the SHG/SFG signal usually enough to not need photon counters¹¹

⁹S. Cattaneo et al. *Optics Letters*, 28(16):144–1447, 2003.

¹⁰F.X. Wang et al. *Physical Review B*, 80(23):233402, 2009.

¹¹P Figliozzi et al. *Physical Review Letters*, 94(4):47401, 2005.

Theory

The nonlinear polarization can be separated into two contributions¹²,

$$\mathbf{P}_{nc}^{(2)} \equiv n_s \left(\gamma^e - \gamma^m - \frac{\gamma^q}{6} \right) (\mathbf{E} \cdot \nabla) \mathbf{E} \equiv |\Gamma_{nc}| e^{i\Phi} (\mathbf{E} \cdot \nabla) \mathbf{E}, \quad (5)$$

$$\mathbf{P}_g^{(2)} \equiv (\delta - \beta - 2\gamma) (\mathbf{E} \cdot \nabla) \mathbf{E} \equiv \Gamma_g (\mathbf{E} \cdot \nabla) \mathbf{E}. \quad (6)$$

The phase Φ causes interference between the glass and nanocrystal signals – this will cause the double peak shape we will see later on.

¹²A. Wirth et al. *Physica Status Solidi C*, 5(8):2662–2666, 2008.

Signal Enhancement

Single beam SHG scales as¹³

$$N_{SHG} \sim \frac{f_{\text{rep}} \mathcal{E}^2}{\tau A^2}, \quad (7)$$

where A is the beam spot size ($A = \pi w_0^2$), τ is the pulse duration, \mathcal{E} is the pulse energy, and f_{rep} is the repetition rate of the pulses.

¹³P. Figliozzi et al. *Physical Review Letters*, 94(4):47401, 2005.

But for two incoming plane wave fields that are cross polarized, SHG counts scale as


$$N_{SHG} \sim \frac{f_{\text{rep}} \mathcal{E}_1 \mathcal{E}_2 \sin^2 \alpha}{\lambda^2 \tau A^2}, \quad (8)$$

where α is the angle between the beams, and λ is the wavelength.

Enhancements

- The $\frac{1}{\lambda}^2$ factor greatly increases signal intensity
- The $\sin^2 \alpha$ term allows us to optimize the beam angle

The $(\mathbf{E} \cdot \nabla) \mathbf{E}$ dependence is directly observable¹⁴:




Experiment (left) and predicted (right).

¹⁴P. Figliozzi et al. *Physical Review Letters*, 94(4):47401, 2005.

Theory

The XP2SHG/SFG Technique

- Three Z-scans needed
- Interference between contributions causes double peak shape¹⁵



¹⁵A. Wirth et al. *Physica Status Solidi C*, 5(8):2662–2666, 2008.

Three huge benefits over single beam

- 1 Intensities are much higher
- 2 Enhanced signal allows better determination of $\chi^{(2)}$
- 3 Dipolar and quadrupolar contributions from nanoparticles can be discerned from substrate bulk contributions

└ The Experiment

└ Laser System, Optics, and Detectors


Laser system output

- Wavelength: 800 nm
- Average power: 1.1 Watts
- Energy: 1 mJ per pulse
- Duration: 100 fs
- Repetition rate: 1 kHz




Noncollinear Optical Parametric Amplifier

- 760 - 515 nm (1.6 - 2.4 eV)
- Energy: 3 - 12 μJ
- Duration: 250 fs
- Repetition rate: 1 kHz



└ The Experiment

└ Laser System, Optics, and Detectors



The NOPA in action.

└ The Experiment

└ Using the XP2SHG/SFG Technique

The XP2SHG/SFG Setup





Diagram of the XP2SHG setup.


└ The Experiment

└ Using the XP2SHG/SFG Technique

XP2SHG using a BBO crystal.




Two beams in. . .



three beams out!


Ellipsometry¹⁶



- Complement nonlinear studies and allow full characterization of nanoparticles
- Ellipsometry can determine material dielectric function and complex index of refraction

¹⁶Image: Stannered

Linear Transmission





Gold samples have plasmon resonance around 530 nm^{17 18}

¹⁷D.M. Schaadt et al. *Applied Physics Letters*, 86:063106, 2005.



¹⁸S. Lin et al. *Advanced Materials*, 17(21):2553–2559, 2005.

Ellipsometry – Substrate



Graphs for Ψ (left) and Δ (right).



Ellipsometry – Gold



Graphs for Ψ (left) and Δ (right). Very noisy, nothing like references¹⁹

¹⁹H.L. Zhang et al. *Advanced Materials*, 15(6):531–534, 2003.

Ellipsometry – Silicon




Graphs for Ψ (left) and Δ (right). Flat transmission curve for Silicon and poor ellipsometry data nothing like previous work.²⁰

²⁰J. Wei et al. *Physical Review B*, 84:165316, Oct 2011.

XP2SHG/SFG – Substrate

Substrate XP2SFG data at $550 + 800$ nm.




Confirms results from previous studies^{21 22}


²¹J. Wei et al. *Physical Review B*, 84:165316, Oct 2011.

²²A. Wirth et al. *Physica Status Solidi C*, 5(8):2662–2666, 2008.

XP2SHG – Gold



Gold XP2SHG data, entry side.




Gold XP2SHG data, exit side.


Characteristics

- Almost identical to substrate data!
- Low input intensity produces no SHG/SFG output
- High input intensity generates white light
- Need to analyze entry and exit positions
- Detector can't distinguish between single and two-beam SHG
- Both dipolar and quadrupolar contribution


XP2SFG – Gold



Noisy signal with no discernible peaks at 520 + 800 nm




White light generation at 560 + 800 nm



Best XP2SFG data at 520 + 800 nm.

- Three different frequencies: NOPA in, 800 nm in, SFG out
- Strong output with NOPA frequency near plasmon resonance
- Wider angle between beams should alleviate scattering

XP2SHG/SFG – Silicon



XP2SFG at 500 + 800 nm.

■ Weaker signal compared to Au

■ Signal not enhanced at all^{a b c}

^aJ. Wei et al. *Phys. Rev. B*, 84, 2011.

^bA. Wirth et al. *Phys. Stat. Sol. C*,
5(8), 2008.

^cP. Figliozzi et al. *Phys. Rev. Lett.*, 94,
2005.

Summary

- Samples of poor optical quality with massive scattering
- Scattering and transparency to blame for erratic linear data
- XP2SHG/SFG technique did not enhance nonlinear signal
- Exciting near resonance may have actually been a hinderance
- Separate Z-scans did not allow for separation of glass and nanoparticle contributions



“So what?”

Conclusions

Accomplished

- Applied the XP2SHG/SFG technique
- Applied complementary linear measurements
- Gained considerable experience

Future Work

- 1 Samples need to be in better shape and better characterized
- 2 Samples mounted on different substrates for linear measurements
- 3 XP2SHG/SFG technique still relatively new with metallic nanoparticles
- 4 Implement this setup here at the CIO

Acknowledgements

- Group leaders, Dr. Ramón Carriles and Dr. Enrique Castro
- Future advisor, Dr. Bernardo Mendoza
- Dr. Mike Downer and Junwei Wei of UT Austin
- Dr. Alejandro Reyes and Dr. Alicia Oliver of the IF-UNAM
- DFA and all CIO staff
- The CONACyT
- Everyone that came to listen to this talk
- Front page image: Travis Jennings
- Second page images: Amanda Barnard, Vulcan 10 Petawatt
- Dr. Cabellos drawn by Luis Adán Martínez

¡Gracias!

¡Feliz cumpleaños a Alberto, Marcelo
&
Edith
y feliz aniversario también!